

## Preparation of Chemically Competent Cells

### TBF I Buffer (pH 5.8)\*

30 mM KOAc  
100 mM RbCl  
50 mM MnCl<sub>2</sub>  
10 mM CaCl<sub>2</sub>  
→ sterile filtered

### TBF II Buffer (pH 7.0)

10 mM MOPS  
10 mM RbCl  
75 mM CaCl<sub>2</sub>  
15 % Glycerol  
→ sterile filtered

## STERILE WORKING TECHNIQUES REQUIRED THROUGHOUT!!!

- 1 colony of *E. coli* into 4 ml LB medium (+ appropriate antibiotics)  
→ 37°C overnight
- 500 µl pre-culture into 200 ml LB medium (+ appropriate antibiotics)  
→ 37°C to OD<sub>600</sub> = 0.45-0.55

## FROM NOW ON WORK ON ICE OR AT 4°C!!!

- collect culture in 4x50mL sterile conical centrifugation tubes (Falcon™ tubes)
- centrifuge at 3000 rpm for 10 min, discard supernatant
- gently (!) resuspend each pellet in 1 ml chilled TBF I and combine in one tube
- fill up to 15mL with TBF I and incubate on ice for 60 min
- centrifuge at 3000 rpm for 10 min, discard supernatant
- resuspend pellet in 4 ml (<sup>1</sup>/<sub>50</sub> of original culture volume) TBF II
- aliquot competent cells (50 or 100 µL) and freeze immediately in liquid nitrogen
- store at -70°C

\* To avoid precipitation of MnO<sub>2</sub> (brown colored solution or precipitate), dissolve everything except MnCl<sub>2</sub> in <sup>3</sup>/<sub>4</sub> of the volume, adjust pH to 6.0 with acetic acid, add MnCl<sub>2</sub> (dissolved in water), fill up with water nearly to the required total volume and then further adjust the pH to 5.8 if required. Avoid adding KOH (or any base) to a solution containing a Mn<sup>2+</sup> salt!!!