

## Quantification of Plasmid DNA

- start the Cary WinUV program “RNA-DNA” and let the UV-lamp warm up
- dilute plasmid DNA from miniprep (1 $\mu$ L DNA + 99 $\mu$ L ddH<sub>2</sub>O) and vortex shortly to mix
- blank UV spec. against ddH<sub>2</sub>O
- measure A<sub>260</sub> and A<sub>280</sub> (make sure that A<sub>260</sub> reading is between 0.1 and 1.0!!!)

$$A_{260} \times 100 \text{ (dilution factor)} \times 50 \text{ ng}/\mu\text{L} = \text{ng DNA}/\mu\text{L} \text{ in undiluted sample}$$

You can expect concentration values between 400 and 500 ng/ $\mu$ L

(Roche plasmid kit, 4mL cell culture, 80 $\mu$ L elution buffer)

### A<sub>260</sub>:A<sub>280</sub> ratio as a purity criteria

--> A<sub>260</sub> Absorption of nucleic acids (and weakly of proteins)

--> A<sub>280</sub> Absorption of proteins (and weakly of nucleic acids)

Rule of thumb: DNA 1.8

RNA 2.0

(Roche kit normally gives you 1.9 for plasmid DNA)