

## Digestion of DNA with Restriction Enzymes, Dephosphorylation of Vector and Purification

### 1) Digestion with NdeI and XmaCI

20  $\mu\text{L}$  vector ( $\sim 200\text{-}500 \text{ ng}/\mu\text{L}$ )  
(OR 30  $\mu\text{L}$  [PCR product](#) after purification on [agarose gel](#))  
5  $\mu\text{L}$  buffer #4 (10x, NEB)  
23  $\mu\text{L}$  (OR 13  $\mu\text{L}$  for PCR product) ddH<sub>2</sub>O, mix well  
2  $\mu\text{L}$  NdeI ( $20 \text{ U}/\mu\text{L}$ , NEB), mix well  

---

50  $\mu\text{L}$

- 18h at 37°C (without shaking)  
(using the PCR machine for this purpose works very well)

- add:

3  $\mu\text{L}$  XmaCI ( $10 \text{ U}/\mu\text{L}$ , NEB)  
0.55  $\mu\text{L}$  BSA (100x, NEB)  
0.5  $\mu\text{L}$  buffer #4 (10x, NEB)  
0.95  $\mu\text{L}$  ddH<sub>2</sub>O, mix well  

---

55  $\mu\text{L}$

- 18-22h at 37°C (without shaking)
- 20min at 65°C (to denature enzymes)

### 2) Dephosphorylation of vector with Shrimp Alkaline Phosphatase (SAP)

55  $\mu\text{L}$  digestion mix  
7  $\mu\text{L}$  SAP buffer (10x, Roche)  
8  $\mu\text{L}$  SAP ( $1 \text{ U}/\mu\text{L}$ , Roche), mix well  

---

70  $\mu\text{L}$

- 2.5h at 37°C (without shaking)
- 20min at 65°C (to denature enzyme)

### 3) Purification of digested vector and PCR product with [agarose gels](#)

- vector: usually  $\sim 1\%$   
PCR products: usually  $\sim 2\%$