


Standard PCR

0.5 μL vector / template	(~200-500 $\text{ng}/\mu\text{L}$)
1 μL forward primer	(100 $\text{pmol}/\mu\text{L}$)
1 μL reversed primer	(100 $\text{pmol}/\mu\text{L}$)
1 μL dNTP mix	(10 mM each)
5 μL 10x polymerase buffer	
41 μL ddH ₂ O, mix well	
0.5 μL polymerase*, mix well	
<hr/>	
50 μL	

- * for longer transcripts and high accuracy use *Pfu* DNA Polymerase (Promega, M774A, 100u, 3u/ μL)
for smaller transcripts or where accuracy is not so important use *Taq* DNA Polymerase (Q·BIOgene, EPTQA025, 250u, 5u/ μL)



- use 0.2mL thin walled PCR tubes
- adjust annealing time to your specific primer if necessary
- adjust elongation (and final extension) time to your transcript length if necessary (rule of thumb: 1:00 for 1000bases, minimum 0:45!)
- for *Taq* DNA Polymerase use an elongation temperature of 72°C