

Transformation of vectors/DNA into chemically competent cells

- thaw 50-100 μL [chemically competent cells](#) on ice
- add 1 μL vector or 10 μL [ligation mix](#) and mix well by gently tipping the tube with your finger
- incubate on ice for 30 min
- perform heat shock in thermomixer (without shaking): 42°C for 45s
(Rosetta and Codon Plus only for 30s!)
- put on ice for 2 min
- add 500 μL pre-warmed (37°C) LB or SOC medium (without antibiotics!)
- shake at $\sim 400\text{rpm}$ for 1h in thermomixer at 37°C
- transformation with vector: plate $\sim 50 \mu\text{L}$ onto an [agar plate](#) containing the appropriate antibiotics
- transformation with ligation mix: collect cells by gentle centrifugation (1000rpm, 2 min), carefully discard supernatant until $\sim 50\text{-}100 \mu\text{L}$ are left in the tube, gently resuspend cells and plate out as above